
# Instructions for use

## Intensiv UniglossPolisher Pat.Pend.EP10162150.6

Scientifically proven and clinically tested at the University of Zurich, Switzerland

### Product description

- Slim flame and compact cup of polymeric material highly loaded with ultra-fine diamonds. The new core shape of the polisher consists of extra resistant plastic and follows exactly the flame and cup shapes.
- Slim flame and compact cup, without colouring.
- Core of flame and cup shape and mandrel consist of one piece of blue-green resistant plastic.
- Right angle instrument, low speed.
- Length of the working part:
  - Flame: 9 mm
  - Cup: 7 mm
- Maximum diameter of the working part:
  - Flame: 3 mm
  - Cup: 6 mm
- Autoclavable, multiple use

### Indications

- Final high gloss polishing of convex and pronounced concave surfaces of aesthetic restorations
- Maintenance of existing aesthetic restorations

### Practical Tips

- Use Intensiv UniglossPolisher on composite restorations after pre-polishing with rotating diamond instruments (15 µm yellow).
- Use water spray during the entire treatment (minimum 50 mL/min).
- Regulate water flow in case of water excess or shortage.
- Apply light force not exceeding 2N.
- Apply speed between 8'000 and 20'000 rpm max.

### Maintenance and Sterilisation

- Instruments are delivered non sterile. They must be sterilized prior to first use and disinfected, cleaned, and sterilized after each use.
- Sterilization must be carried out according to validated procedures.
- For further information please refer to "Maintenance and Sterilization guidelines" chapter.

### Risk Warnings

- Temperatures above 140°C should be avoided as they could increase the risk of mandrel fracture during usage.
- Avoid operating force of over 2N, because of fracture risk of the plastic shaft and of the tip in the flame shaped polisher.
- To ensure traceability of the instruments during their entire application, we recommend keeping the packaging.

### Maintenance and Sterilization guidelines

- Side effects: None to our knowledge.
- Working place: Remove all dirt on surfaces using a paper towel.
- Storage and transportation: No special requirements.
- Pre-cleaning steps: See "Disinfection" and "Sterilization". No dismantling needed.
- Automatic cleaning: Material required: cleaning/disinfecting machine and cleaning product.
  - 1) Place instruments in the machine.
  - 2) Set the machine on standard cycle and clean for at least 5 minutes and rinse for at least 3 minutes.
- Disinfection: A disinfection solution can be used, according to instructions on label (Dürr Dental ID 212, for instance). In case of automatic cleaning, 10 minutes rinsing at 80°C can be performed as last step.
- Maintenance: Discard any blunt or damaged polisher.
- Checking of the instruments: For all instruments: ensure that the instrument is neither damaged nor worn. Its' edges and corners must have no notches and be regular in shape.
- Sterilization: Use a single-pulsed or fractionated vacuum autoclave and subvacuum drying. Temperature must not exceed 140°C. Do not use a chemiclave or a hot air disinfecter.
- Additional information: If several instruments are to be sterilized in the autoclave, ensure that the maximum weight allowed by the autoclave manufacturer is not exceeded.

Intensiv UniglossPolisher,  
trial pack with 3 polishers each, flame and cup  
Ref. 2300/6

Intensiv UniglossPolisher, pack of 6 polishers, flame  
Ref. 2310/6

Intensiv UniglossPolisher, pack of 6 polishers, cup  
Ref. 2320/6

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