

# Intensiv ProxoStrip Pat. pend. 2007 0219/07

A unique diamond-coated metal strip  
for the final polishing of proximal surfaces  
**Practical, rapid, ergonomic!**

- Finishing and polishing in one step
- Ergonomic grip for a simple and easy use
- Short strip for an improved finger support while stripping
- Clinically perfectly polished proximal surfaces
- Sterilizable, reusable

Intensiv ProxoStrip is a unique diamond-coated metal strip with two working parts in different grit sizes (40µm red and 15µm yellow) and perforated ends for a safe grip. The strip height of 2,5mm as well as the uncoated zone between the two working parts preserve the contact point.

## Indications

- Finishing and polishing of fillings and crown borders in proximal areas
- Removal of excess material and overhangs in fillings and crown borders
- Bevelling in adhesive restorations in cervical proximal areas



For the bevelling of preparation margins both the red and the yellow working part can be used.



ProxoStrip is also suitable for a proximal finishing without damaging the contact point.



The yellow working part is well-suited for a clinically perfectly polished proximal surface without impairing the contact point.

Ref. PX 4015-6  
Package of 12 Strips



Ref. PX 4015-12  
Package of 12 Strips



Ref. PX 4015
µm
■ 40
■ 15

Intensiv SA  
via al Molino 107  
P.O. Box 741  
6916 Grancia  
Switzerland

Tel. + 41 91 986 50 50  
Fax + 41 91 986 50 59  
info@intensiv.ch  
www.intensiv.ch

