

Intensiv HygienicTray

by M. Gigandet, T. Engel, University of Bern

The unique bur block developed
by dentists for dentists!

Hygienic, compact, sterilizable!

Bur blocks have become indispensable for the proper organisation of rotary instruments in any well-organized dental practice.


- Easy to clean due to the complete dismantling in 3 parts (cover, instrument tray, bottom)
- Ergonomic compact design, suitable for any standard sterilisation tray
- Ideal arrangement of 13 FG and 5 RA holes for easy access to the instruments
- Resistant, high-quality stainless steel

Indications

- Proper organisation of the instruments during the treatment
- Support for a specific indication treatment (e.g. Restorative bur block)
- Bur holder during the sterilisation process



Ref. HT 100



Intensiv SA
via al Molino 107
P.O. Box 741
6916 Grancia
Switzerland

Tel. + 41 91 986 50 50
Fax + 41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

 **Intensiv**
Swiss Dental Products