

Intensiv ApproxOpener

Manual and serrated metal strip for removal of coronal excess of adhesive aesthetic restorations



Indications

- Removal of emerged adhesive material in the coronal area

Benefits

- Removal of emerged adhesive material without scratches thanks to not diamond-coated version
- Secure strip movement thanks to tension of the strip in the bow
- Sterilisable, reusable

Product Description

- Manual metal strip clamped into a small metal bow
- Serrated
- Height of strip: 4 mm
- Length of strip: 20 mm



REF AO2018/3