

Intensiv Cerinlay Set

Special diamond instruments for inlay cavity preparations

Restorations using ceramic inlays require simple preparations with rounded interior edges and sharp cavity margins without bevelling. The slightly tapered shape and suitable small diameters of these instruments make undercut-free cavity preparations possible.



REF 011

Product Description

- 4 cones, rounded edges, Ø 014 and 018, 6.0 and 8.0 mm, 80 µm
- 4 shape and dimension matching instruments in 25 µm

Colour code	Micron	Description
■ 524	80 µm	Medium
□ 514	25 µm	Fine

REF 011 - FG	Instrument	R	L	Ø	Colour Code
■	FG 8425	0.28	6.0	014	Medium
	FG 8427	0.20	8.0	014	Medium
■	FG 8525	0.28	6.0	018	Medium
	FG 8526	0.28	8.0	018	Medium
□	FG 3425	0.28	6.0	014	Fine
	FG 3427	0.20	8.0	014	Fine
□	FG 3525	0.28	6.0	018	Fine
	FG 3526	0.28	8.0	018	Fine

Indications

- Preparation and finishing of ceramic inlay cavities for premolars and molars

University of Berlin, Germany



Benefits

- Preparation method clinically tested at the University of Berlin
- Finishing of fracture-free enamel margins
- Split-free marginal integration of inlays

Clinical pictures

- 01 Ceramic inlay cavity preparation with instrument FG 8526
- 02 Cavity finishing with instrument FG 3526
- 03 Cavity finishing with instrument FG 3525