

Intensiv Ortho-Strips® Sets

Neu

Set für approximale Schmelzreduktion

Das Intensiv Ortho-Strips Tray 400 Set ist für die ASR (Approximale Schmelzreduktion). Es enthält auch die neuen, patentierten Intensiv Ortho-Strips System, Central. Dank der Diamantbeschichtung im zentralen Arbeitsbereich, bietet das Intensiv Ortho-Strips System, Central eine hohe Präzision bei der Reduktion sowie dem approximalen Finieren und Polieren des Zahnschmelzes unter Respektierung der ursprünglichen Morphologie und des Zahngewebes.



REF OST400Set01/6



REF OST400Set02/6



REF OST400Set03/6

Produktbeschreibung

REF OST400Set01/6

- 1 Intensiv Ortho-Strips Tray REF OST400
- 1 Intensiv Ortho-Strips System, Opener, REF OS08OP-DS
- 2 Intensiv Ortho-Strips System, Central, Medium, REF OS40M-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Coarse, REF OS60C-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Fine, REF OS25F-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Polishing, REF OS15POL-CEN-DS

ERWEITERT

Öffnen (Orange), konturiertes reduzieren (Rot) und rasches Reduzieren (Braun), Finieren (Weiss), Polieren (Gelb)

REF OST400Set02/6

- 1 Intensiv Ortho-Strips Tray REF OST400
- 1 Intensiv Ortho-Strips System, Opener, REF OS08OP-DS
- 1 Intensiv Ortho-Strips System, ExtraCoarse, REF OS80XC-DS
- 1 Intensiv Ortho-Strips System, Central, Coarse, REF OS60C-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Medium, REF OS40M-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Fine, REF OS25F-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Polishing, REF OS15POL-CEN-DS

KOMPLETT

Öffnen (Orange), konturiertes reduzieren (Rot) und rasches Reduzieren (Braun), sehr rasches Reduzieren (Blau), Finieren (Weiss), Polieren (Gelb)

REF OST400Set03/6

- 1 Intensiv Ortho-Strips Tray REF OST400
- 2 Intensiv Ortho-Strips System, Opener, REF OS08OP-DS
- 3 Intensiv Ortho-Strips System, Central, Medium, REF OS40M-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Polishing, REF OS15POL-CEN-DS

STANDARD

Öffnen (Orange), konturiertes reduzieren (Rot), Polieren (Gelb)

Indikationen

- Set03: Ortho-Strips Standard-Protokoll; Öffnen, Reduzieren/Messen und Polieren
- Set02: Ortho-Strips Vollständige Auswahl
- Set01: Ortho-Strips Standard-Protokoll inkl. 60µm (braun) für raschere Abrasion und 25µm (weiss) für anschließende Finitur

Vorteile

- Die Intensiv Ortho-Strips sind während der Behandlung gemäss Zuordnung ihrer Anwendung verfügbar
- Vereinfachtes Einführen durch den interdentalen Kontaktpunkt
- Vermeidung von zervikalen Stufenbildungen
- Dentinüberempfindlichkeiten werden nicht gefördert