

Intensiv Burs-SizeCard

New!

Unique gauge card with continuous scale

Practical stainless steel card to measure dimensions of rotating instruments

Dental clinics have the necessity to measure or to control the lengths and diameter of rotating instruments in use.

Indications

- Diameter validation of rotating instruments during treatment
- Identification of instrument dimensions in case of package label lack

■ Length scale for measuring the instrument total length and working part: from 0 to 70mm

Intensiv Burs-SizeCard
Ref. 6100




■ Continuous slots for diameter measurement of rotating instrument

■ Diameter measurement scale: 006 to 038 (1/10mm)

1:1

- Precise definition of rotating instrument diameter and length
- Proper re-ordering of the instrument used in the dental clinic
- Measurement of rotating instrument exact dimension (+/- 0.01 Tolerance) due to the continuous scale



Intensiv SA
6926 Montagnola
Switzerland

Tel. +41 91 986 50 50
Fax +41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

 **Intensiv**
Swiss Dental Products