

New!

Intensiv HygienicTray HT3000

Developed by dentists for dentists

Resistant high quality stainless steel instrument tray for the disposal of all rotating instruments for routine treatment protocols

Intensiv HygienicTray HT100

- Folded back cover serves as base for the tray during treatment
- Ergonomic instrument access due to high and inclined tray position



Ref. HT100*

Intensiv HygienicTray HT200


- Ergonomic instrument access to stand-alone tray



Ref. HT200*

Intensiv HygienicTray HT300

- Folded back cover keeps the instrument tray in place during treatment
- Ergonomic instrument access due to low and inclined tray position



Ref. HT300*

Indications

- Organization of standard instruments for all routine treatment protocols
- Instrument holder during the sterilization process
- Storage of rotary instruments

Benefits

- Organization and availability of commonly used rotating instruments
- Easy instrument removal due to irregular arranged position of instruments
- Dismountable into 2 parts – free from fissures and hollows – completely hygienic
- Durable resistant without limitation sterilizable stainless steel



Ref. HT3000*

*delivered without rotating instruments

Intensiv SA
6926 Montagnola
Switzerland

Tel. +41 91 986 50 50
Fax +41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

 **Intensiv**
Swiss Dental Products