

Intensiv ProxoContour^{Pat. 699 819}

Clinically tested by the University of Zurich, Switzerland

New!

The only in two grit sizes one sided diamond-coated metal strip with ergonomic handles for contouring and finishing of proximal aesthetic restoration surfaces

Benefits

- Efficient manual contouring of the proximal filling morphology
- Finishing of the roughened surfaces in one process
- Short strip for improved finger support during treatment
- No injury risk of soft tissue (lip cutting)
- Sterilizable, reusable

Intensiv ProxoContour metal strip with one sided diamond-coated grit sizes: 60 µm (brown) for contouring and 40 µm (red) for finishing.

The 2.5 mm height of the strip allows the secure stripping below the contact point.

The diamond-free zone between the diamond coated surfaces provides a non-destructive contact zone passage.

Indications

- Contouring and finishing of aesthetic tooth fillings in the proximal area
- Removal of excess material and contour on proximal restoration surfaces



Clinical pictures:
Dr. Simona Giani, Varese, Italy



Pack of 6 strips
Ref. PXC6040/6

