

Intensiv ProxoPolish

Pat. 699 819

Clinically tested by the University of Zurich, Switzerland

New!

The only in two grit sizes one sided diamond-coated metal strip with ergonomic handles for efficient stain removal and manual polishing of proximal surfaces

Benefits

- Efficient manual cleaning and polishing in one treatment process with the same instrument
- Ergonomic strip handles and short strip for simple and easy use
- No injury risk of soft tissue (lip cutting)
- Sterilizable, reusable

Intensiv ProxoPolish is the only one sided diamond-coated metal strip in two grit sizes 15 μm (yellow) and 8 μm (orange).

The 2.5 mm height of the strip allows secure stripping below the contact point.

The diamond-free zone between the diamond coated surfaces provides a non-destructive contact point passage.


Indications

- Stain removal in proximal areas
- Gloss polishing in proximal areas
- To complete the professional dental hygiene in the proximal area

Clinical pictures:
Dr. Simona Giani, Varese, Italy



Pack of 6 strips
Ref. PXPO1508/6



Intensiv SA
6926 Montagnola
Switzerland

Tel. +41 91 986 50 50
Fax +41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

Intensiv
Swiss Dental Products