

Intensiv ProxoPolish

Pat. 699 819

Clinically tested by the University of Zurich, Switzerland

New!

The only in two grit sizes one sided diamond-coated metal strip with ergonomic handles for efficient stain removal and manual polishing of proximal surfaces

Benefits

- Efficient manual cleaning and polishing in one treatment process with the same instrument
- Ergonomic strip handles and short strip for simple and easy use
- No injury risk of soft tissue (lip cutting)
- Sterilizable, reusable

Intensiv ProxoPolish is the only one sided diamond-coated metal strip in two grit sizes 15 μm (yellow) and 8 μm (orange).

The 2.5 mm height of the strip allows secure stripping below the contact point.

The diamond-free zone between the diamond coated surfaces provides a non-destructive contact point passage.


Indications

- Stain removal in proximal areas
- Gloss polishing in proximal areas
- To complete the professional dental hygiene in the proximal area

Clinical pictures:
Dr. Simona Giani, Varese, Italy



Pack of 6 strips
Ref. PXPO1508/6



Intensiv SA
6926 Montagnola
Switzerland

Tel. +41 91 986 50 50
Fax +41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

 **Intensiv**
Swiss Dental Products

Intensiv Innovations


Please send your order via fax to 0041 (0) 91 986 50 59 or via mail to order@intensiv.ch

Product	Ref.	Price €	Quantity
Intensiv UniglossCellbrush	2200/6	28,80	
Intensiv UniglossCellbrush	2200/12	54,00	
Intensiv UniglossPaste Syringe	2500	38,70	
Intensiv UniglossPaste Blister 25	2510	89,00	
Intensiv ProxoStrip	PX4015/6	44,00	
Intensiv ProxoStrip	PX4015/12	70,55	
Intensiv ProxoStrip Plus	PXP4015/6	49,00	
Intensiv ProxoStrip Plus	PXP4015/12	78,40	
Intensiv ProxoContour	PXC6040/6	55,20	
Intensiv ProxoPolish	PXPO1508/6	57,00	
Intensiv HygienicTray	HT100	46,80	
Intensiv HygienicTray	HT200	46,80	
Intensiv HygienicTray	HT300	33,60	
Intensiv HygienicTray	HT3000	90,00	
Intensiv PerioDiaCurette	PDC-U	53,90	
Intensiv PerioDiaCurette	PDC-11/12	53,90	
Intensiv PerioDiaCurette	PDC-13/14	53,90	
Intensiv PerioDiaCurette	PDC-F	53,90	
Intensiv PerioDiaCurette	PDC-A	202,70	
Intensiv MarginFinisher	MF40A	50,80	
Intensiv MarginFinisher	MF40P	50,80	
Intensiv MarginFinisher	MF40A+P	95,05	
Intensiv Proxoshape Flexible	PSF15	24,60	
Intensiv Proxoshape Flexible	PSF15L	24,60	
Intensiv Proxoshape Flexible	PSF15/6	140,20	
Intensiv Proxoshape Flexible	PSF15L/6	140,20	
Intensiv Proxoshape Flexible	PSF40	24,60	
Intensiv Proxoshape Flexible	PSF40L	24,60	
Intensiv Proxoshape Flexible	PSF40/6	140,20	
Intensiv Proxoshape Flexible	PSF40L/6	140,20	
Intensiv Proxoshape Flexible	PSF60	24,60	
Intensiv Proxoshape Flexible	PSF60L	24,60	
Intensiv Proxoshape Flexible	PSF60/6	140,20	
Intensiv Proxoshape Flexible	PSF60L/6	140,20	

My dental clinic:

Date _____ Signature _____

Dental Dealer:



Please send an Intensiv catalogue

Intensiv SA
6926 Montagnola
Switzerland

Tel. +41 91 986 50 50
Fax +41 91 986 50 59
info@intensiv.ch
www.intensiv.ch

