

Intensiv ProxoContour Coarse

Pat. 699 819

Clinically tested by the University of Zurich, Switzerland

New!

Diamond coated and perforated metal strip with grips, for modeling and finishing

During creation of pronounced direct and indirect adhesive restorations undesirable material excess on approximal surfaces and margins could occur. Material excess must necessarily be removed to ensure the stability of the restoration itself and its long term clinical success.

Indications

- Removal of marginal material excess of direct and indirect esthetic reconstructions
- Contouring and finishing of specific excess on the approximal surfaces of pronounced aesthetic adhesive restorations

Benefits

- Friction-free application due to perforation of the diamond coated strip part
- Fast and efficient removal of excess material due to the coarse grain size
- 2 grips for improved manual handling during treatment



Package of 6 strips
Ref. PXCC8040/6

To sterilize, multi-use



Clinical pictures:
Dr. A. Devigus, Bülach, Switzerland

